
# Operating Instructions

Household  
Electric Toothbrush

**EW118**

**Panasonic**

Before operating this unit, please read these instructions completely.

# CAUTIONS REGARDING SAFETY

- Before using this product, please read all cautions thoroughly and use the product properly
- Following the cautions stated herein will ensure that you use the product safely and will prevent harm or injury to yourself and others
- The harm or injury which may occur when the product is not used properly is also described in the cautions

The cautions contain important information regarding safety and must be followed

Improper use may result in personal injury or product damage as described below

## CAUTIONS

- Do not press the brush against the gums too strongly. Doing so may result in damage to the gums. Persons with severe pyorrhea should consult their dentist
- Do not discard batteries in fire. Doing so may cause the batteries to explode, resulting in injury or damage
- Do not attempt to repair or modify the product. Doing so may cause the product to malfunction, resulting in injury
- Close supervision is necessary when this product is used by, on or near children or invalids
- Use this product only for its intended use described in this manual. Do not use attachments not recommended by the manufacturer
- People with tooth or gum disease and people receiving oral treatment should consult a dentist before use
- Do not open the power handle except to insert or remove battery. Damage to the product and/or personal injury may occur
- **KEEP FLOSS ATTACHMENTS AWAY FROM SMALL CHILDREN**. They are small enough to be swallowed and can become lodged in the throat
- Never use while drowsy

Please save these cautions for later reference

fig. ①


fig. ②





fig. ③

**fig. ④**

Horizontal Bass/Floss Motion  
(when held horizontally)

Use side edge at gumline

Use heel edge


Use toe edge



\* e

**fig. ⑤**

Horizontal Bass/Motion  
(when held horizontally)

c (clean between teeth)


d (slowly and gently up & down)

SIDE VIEW

e (curve floss around tooth in C-shape)


TOP VIEW


\* Drawings taken from

Wilkens, E M "Clinical Practice of the Dental Hygienist" 6th Edition Philadelphia, Lea & Febiger, 1989, pp 304 - 305

**fig. ⑥**

**fig. ⑦**

## Parts Identification (See fig ①)

- Ⓐ Power handle
- Ⓑ Cap (Removable when cleaning)
- Ⓒ Metal shaft (Insert toothbrush onto shaft)
- Ⓓ On/off switch
- Ⓔ Battery cover
- Ⓕ Toothbrush (1)  
(Small size included)
- Ⓖ Flosser accessories
- Ⓗ Flosser stem (1)
- Ⓘ Floss attachments (30pcs)
- Ⓢ Storage case

## Before use

### Battery insertion (See fig ②)

#### NOTE

Use **ONLY AA (LR6) alkaline battery**

- The power handle cannot operate at full power if you use other types of batteries
  - When inserting or replacing the battery, wipe dry power handle to ensure that no water get inside
- 1 Turn the battery cover counterclockwise to remove
  - 2 Insert one **AA (LR6) alkaline battery** in the power handle as shown. Then re-attach the battery cover firmly, turning in clockwise direction

## How to use

### Using your toothbrush (See fig ③)

- 1 Attach toothbrush by inserting hollow end of brush into metal shaft. Press down until it snaps firmly in place.
- 2 Move the on/off switch to "on". With the new electric toothbrush, undue force should not be applied to teeth or gums.
- 3 When using toothpaste, put brush head into your mouth before turning on switch to avoid splattering of toothpaste.

### After brushing

- Rinse bristles with water and shake drip-dry
- Rinse brush handle with water and wipe dry
- Wash brush once weekly
- Never submerge handle in water

## Maintenance

- 1 After each use, detach the toothbrush from metal shaft and rinse under running water. **DO NOT USE HOT WATER** (See fig ⑥). Shake off excess water, then store toothbrush and power handle in storage case.
- 2 To remove excess toothpaste and other debris inside cap and around rubber seal, squeeze notched side of the cap (B) and pull upward to remove. Do not nick waterproof rubber part. Clean with soft cloth after running under water (See fig ⑦). Upon completion, attach cap all the way down. Power handle may be completely immersed in water for quick cleaning. Dry thoroughly after immersion, prior to use.
- 3 To clean power handle and storage case, use damp soft cloth with mild soapy water. Do not wipe with benzene, paint thinner or any other solvents.

### Horizontal (Bass) Motion (See fig ④)

This brushing motion is effective in removing plaque under or along gumline

- 1 Place bristle on gumline at a 45° angle and brush for 3 ~ 5 seconds at each point as shown in ④ a
- 2 Use heel edge for inside surfaces of front teeth as shown in ④ c
- 3 Use toe edge for the very back side of backteeth (last molar) as shown in ④ d
- 4 Apply gentle pressure since excess pressure may injure your gum. If you suffer from severe pyorrhea, consult your dentist before use
- 5 When moving the brush to an adjacent position, overlap the brush position as shown in ④ e
- 6 As a guide, try to brush 2 to 3 minutes, covering all areas & utilizing both motions for most efficient brushing

#### Pointers

- 1 Any type of toothpaste can be used but excessive amount should be avoided if possible
- 2 Using mirror while brushing may help
- 3 If bristles are bent outward, throw out and replace with new brush

### Power Floss (See fig ⑤)

Using power floss is effective between teeth, especially where brush cannot reach

- 1 Attach floss attachment to flosser stem by sliding in. Make sure attachment "clicks" in ⑤ a
- 2 Insert flosser assembly to metal shaft on power unit ⑤ b
- 3 Move the on/off switch to activate movement
- 4 Slowly and gently follow contour of inner side of each tooth surface as shown ⑤ c d e f & g. It may be necessary to move floss up & down a couple of times to clean plaque off tooth. DO NOT press down hard. It may damage or cut gum
- 5 Repeat above process for all inner areas between teeth
- 6 Rinse floss in water when visible debris is on floss, and/or wipe off with clean tissue paper
- 7 After flossing all areas, clean floss stem, attachment and power handle thoroughly by rinsing under running water. DO NOT USE HOT WATER. Remove excess water droplets before reinserting into storage case
- 8 Disposal of attachment after each use or once floss string becomes frail is recommended
- 9 For hygienic reasons different people should never use the same floss string. Please replace it each time an individual uses it

## **WARNINGS (Regarding Floss)**

- a Person with mild to severe periodontal disease SHOULD CONSULT WITH DENTIST before trying the device
- b Children under 3 years old MUST NOT USE this device
- c Children between the ages of 3 to 7 must not use device unless parent does the brushing and flossing for them USE DEVICE ON CHILD ONLY AFTER CONSULTING WITH DENTIST
- d Children over 7 years old must use under strict supervision of parent ONLY AFTER CONSULTING WITH FAMILY DENTIST
- DO NOT press down hard on gums when flossing ALWAYS READ DIRECTIONS BEFORE USING DEVICE
- Some bleeding may occur when initially using the flossing device Do not be alarmed If bleeding persists for more than 7 consecutive days, stop using floss and consult with your dentist

## **NOTE**

- If for any reason floss attachment comes off the flosser stem and gets caught in between teeth, push middle STOP button immediately Remove floss attachment slowly from mouth ONLY after motion stops completely
- If floss string seems to break at same spot, stop flossing that section It may mean that there is an obstruction such as a filling that the string gets caught on (See your dentist to remedy the problem)
- Make sure main body is held perpendicular to tooth surface as shown (If floss is angled incorrectly, it may be difficult to insert or remove as shown in ⑤ h)
- Never use flosser as shown in ⑤ i Flosser MUST NOT cross over a row of teeth Always use as shown ⑤ j
- Using mirror while flossing may help

## Specifications

<b>Power source</b>	DC 1.5V (AA alkaline battery, not included)
<b>Speed</b>	Approx. 3,000 rpm (no load)
<b>Weight</b>	Power handle 3.5 Oz (100g)
<b>Dimensions</b>	Power handle 15/16" (Diameter) x 6-47/64" (Length) 24mm (Diameter) x 171mm (Length)
	Storage case 1-3/16" (H) x 6" (W) x 2-1/16" (D) 30 (H) x 178 (W) x 53 (D)mm
<b>Accessories:</b>	Toothbrush Small size x 1 pc
	Floss Floss attachments x 30 pcs
	Flosser stems x 1 pc
	Storage case

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